

TRAPEZE BAR



The Trapeze Bar is designed to assist the patient with a means of self-help to change positions in bed, to move onto a bedpan, to move from the bed to a bedside commode, or to transfer to and from a wheelchair with minimal help from an attendant.

Patients and caregivers should give careful consideration to the adjustment of the grab bar position before the technician leaves your home. This will prevent the necessity for you to loosen the clamps to readjust the trapeze bar later.

All clamps should be checked daily to ensure that they are securely tightened.

When the head of the bed is elevated or when a caregiver is working with the patient in bed, the triangular grab bar can be placed out of the way by hooking it over the upper end of the off-set bar and sliding it back towards the head of the bed.

IMPORTANT: Please discontinue use and call our office immediately at 800-473-4669 in the event the Trapeze Bar ever becomes loose or unstable in any way.